



E. coli
Detected!

TECTA™ - When Every Second Counts Why Wait For Results?

When regularly testing water quality you want to be certain you get the results back quickly, accurately and efficiently.

Fast Accurate Results in 2-18 Hours

The TECTA™ B16 is fully automated and can provide *E. coli* and Total Coliform results in 2-18 hours, depending on the level of contamination. With integrated networking capabilities, the TECTA™ B16 will provide immediate notification and early warning of positive sample results as soon as they occur. Results are automatically sent to any device, including computers, tablets and smart phones, upon detection of bacteria.

***E. coli* Detected - When the health of many is at stake, you need answers fast!**

